

Standing on The Shoulders of Giants: Metabolomics Made Easy for Non-Experts


By Adopting Community Best Practices with Intelligent Automation

J. Will Thompson, James Campbell, J. Scott Mellors, and J. Michael Ramsey
Move Analytical LLC, Carrboro NC

Introduction

Democratizing access to quantitative metabolite measurements is key to expanding metabolomics into new fields and gaining new users. Currently, most non-experts access metabolomics through service models due to steep learning curves across sample prep, method development, data collection, and analysis. This poster presents an end-to-end workflow using CE-MS, consumables, and intuitive software—all aligned with community best practices—to make hands-on metabolomics accessible to non-experts.

The pillars of an accessible metabolomics workflow are outlined below. The workflow is being sold as MoveKit™ CE.



These pillars can be expanded to other LC-MS based workflows.

Consumables for Rapid Method Deployment



What's Included:
Extraction and collection 96-well plates
Extraction and Focusing Reagent
Internal standards (36 SIL)
SST, Blank, Calibration, and QCs
Software

Protocol
20 µL sample (10:1 diluted for media)
180 uL of extraction+focusing reagent
Shake (15 min)
Filter (spin or vacuum manifold)

Coverage:
Library of 300+ polar metabolites
14 metabolite classes; no derivatization
Flexibility to add your own compounds to the method!


Reduce Decision and Training Fatigue.
Automate the Pain Points.
Unlock Data Exploration.

Materials and Methods


Metabolites were extracted from 20 µL sample using the MoveKit™ CE kit (Move Analytical), which includes reagents, plates, stable-isotope internal standards (Cambridge Isotope Labs), calibration standards, and QCs. The protocol (protein precipitation + filtration) enables 96-well processing in under an hour. Samples were analyzed using Peptides BGE and HR chip on the ZipChip CE system (Repligen) coupled to an Exploris 240 MS (ThermoFisher). The MoveDesk software guided users through experimental setup, QC, and data visualization. Skyline enabled transparent data analysis and quantification. The MoveDesk™ software package was written using Electron, React, Python, and FastAPI.

Software Workflow and Automation


An intuitive UI allows users to make critical experimental design decisions quickly, review the acquisition queue and plate layout, and produce instrument-ready analysis sequences for the batch (below).



System Suitability Test (SST) results are produced automatically. During the data collection, the Data Check feature produces real-time feedback on data quality including intensity, migration index, and separation resolution.



The Analysis uses a Skyline template with automated peak assignment correction. Following quantification, filtering on precision of Study Pool Quality Controls (SPQC), % missing data, and PCA for outlier detection allows only the highest quality measurements to be passed along for statistical comparison or visualization during Reporting.



Acknowledgements and Conflicts

Move Analytical is a startup company focused on developing tools to **Make Metabolomics Easy**. All authors are cofounders and own interest in Move Analytical. We gratefully acknowledge our continued collaborations with Dr. Andrew Percy (Cambridge Isotopes), Dr. Whitney Stutts (UNC), Dr. Laura Herring (UNC), Dr. Sean Mack (NIH), Dr. Jeremy Wolff (NIH), Dr. Erik Allman (AstraZeneca), Dr. Cristina DiPoto (AstraZeneca), and Dr. Erin Redman (Repligen Corp) in the development of our platforms.



www.moveanalytical.com
info@moveanalytical.com

For Research Use Only, Not for Diagnostic Use.
The content of this poster is the subject of one or more patents pending.